-  *Bedienungsanleitung*
-  *Instruction Manual*
-  *Manual de instrucciones*
-  *Mode d'emploi*

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
ECO-MATIC®

by HSW

➤ *Automatic self-filling syringe for oral and pour-on application or injection*

- A) Ergonomically designed handle
- B) Plastic barrel, piston, O-ring and felt lubrication ring
- C) Fixing screw
- D) Drench-nozzle
- E) Hanging Hook
- F) Dosage Adjustment Wheel
- G) Tubing attachment


Instruction manual

The 12.5ml, 30ml and 70ml HSW ECO-MATIC® were developed for oral and pour-on application of solutions for sheep and cattle or multiple dosing of live stock. The ergonomic and user-friendly design, as well as production with high strength materials ensures fireless use as well as a long service life.


It is essential that the syringe is cleaned before each use. Please follow the manufacturers' recommendations with any cleaning or disinfectant solutions used.

Nozzle sizes for HSW ECO-MATIC® 12,5ml, 30ml

All accessories for drenching, pour-on or injection provided for the HSW ECO-MATIC® 12.5 ml and 30 ml model have the same size. Thus, the different nozzles are easily interchangeable between these devices.

Nozzle sizes for HSW ECO-MATIC® 70ml


Accessories provided for the HSW ECO-MATIC® model with 70 ml volume show a larger diameter at the nozzle and the outlet of the syringe. This implies that accessories delivered for HSW ECO-MATIC® 70 ml models are not interchangeable with accessories for the HSW ECO-MATIC® 12.5 ml and 30 ml model.

1) To adjust the dosing amount

Turn the dosing wheel clockwise or counter-clockwise until the desired dose is indicated in the window.

2) Preparing the HSW ECO-MATIC® for use

1. Fix tube on tubing attachment (G) of the HSW ECO-MATIC®.
2. Option 1: Fix the tube on the attachment of the plastic backpack.
Option 2: Fix the tube on the draw-off cap for use with bottle.
3. The drench-nozzle can be adjusted to any angle based on the application by loosening the fixing screw and rotating the nozzle.
4. Adjust HSW ECO-MATIC® to the maximum dosage with the nozzle pointed upward at a 45 degree angle.
5. Depress and release handle repeatedly to draw solution into the plastic barrel.
6. If there are air bubbles, repeat the procedure until there are no more.
7. Please pay attention to the fact that the HSW ECO-MATIC® is completely filled and that there are no air bubbles in the plastic barrel. If necessary, turn the drench-nozzle upwards so that air can escape.
8. The device is now primed and can be adjusted to any dosage as necessary.

As accessories, further metal and plastic nozzles as well as Luer Lock adapters for the use as an injection syringe, can be delivered.

3) Use as an injection syringe (only HSW ECO-MATIC® 12,5 ml and 30 ml)

For this application special accessory is necessary:

1. Remove cap nut and drench nozzle.
2. Insert Luer Lock attachment and replace cap nut.
3. Insert Luer Lock needle.



Attention:

- Leave the protective cover (only on disposable needles) on the needle until ready for use.
- Insert needle onto Luer Lock fitting and secure the needle clockwise.
- Disconnect the needle fitting the opposite way. Replace protective cover for disposal.

INSTRUCTION AFTER USE OF THE HSW ECO-MATIC®

Remove all traces of material from syringe by flushing thoroughly with a hot water detergent mix. Rinse thoroughly with clean, warm water.


Change/cleaning of the valve

- Remove the cap nut from the cylinder by rotating
- Take off the accessory
- fig.1: For HSW ECO-MATIC® with 12,5 ml and 30 ml volume: take-off the valve with spring and O-ring ahead against a little resistance. Do not pull at the spring!
- fig.2: For HSW CO-MATIC® with 70 ml volume, take out the valve with spring and O-ring by screwing out the valve screw at the front end of the barrel.
- Change or clean the valve with spring and O-ring.
- fig. 1: For HSW ECO-MATIC® with 12.5m and 30 ml volume, put the valve incl. spring and O-ring into the barrel and let it click into place against a little resistance.
- fig. 2: For HSW ECO-MATIC® with 70 ml volume, put the valve incl. spring and O-ring from the front into the barrel and fix it by screwing the valve safety screw.

Please note that the valves, springs and O-rings of the HSW ECO-MATIC® 12.5 ml and 30 ml are not interchangeable with valves, springs and O-rings of the HSW ECO-MATIC® 70 ml.

Cleaning the valve of the pour-on Nozzle for the HSW ECO-MATIC® 70ml

The outlet valve of the pour on nozzle for the HSW ECO-MATIC® 70 ml model can be cleaned after unscrewing the spray nozzle at the front. For proper reassembling the valve has to be inserted as shown on the picture on the left.

DISINFECTION

It is recommended to disinfect the HSW ECO-MATIC® and accessories before and after each use.

- Pump hot water through the HSW ECO-MATIC®.
- Suspend the HSW ECO-MATIC® and accessories used for application of solutions in a container of water and boil for at least 20 minutes
- Remove the HSW ECO-MATIC® from container, wrap a cloth around the handle and pump dry.

TO LUBRICATE:

Remove plastic barrel by rotating. Lubricate O-ring and felt ring with 100% pure olive oil, cold squeezed or equivalent vegetable oil. Please follow the manufacturer's durability and storage life instruction. Do not use mineral oil, silicone, WD-40, glycerine or alcohol. Petroleum based products harden the O-rings, reducing the life of the component.

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